

---

# Microbiology of food and animal feeding stuffs — Horizontal method for the enumeration of sulfite- reducing bacteria growing under anaerobic conditions

ICS 07.100.30











## Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO 15213 was prepared by Technical Committee ISO/TC 34, *Food products*, Subcommittee SC 9, *Microbiology*.

## Introduction

Because of the large variety of food and feed products, this horizontal method may not be appropriate in every detail for certain products. In this case, different methods specific to these products may be used if absolutely necessary for justified technical reasons. Nevertheless, every attempt will be made to apply this horizontal method as far as possible.

When this International Standard is next reviewed, account will be taken of all information then available regarding the extent to which this horizontal method has been followed and the reasons for deviations from this method in the case of particular products.

The harmonization of test methods cannot be immediate, and for certain group of products International Standards and/or national standards may already exist that do not comply with this horizontal method. It is hoped that when such standards are reviewed they will be changed to comply with this International Standard so that eventually the only remaining departures from this horizontal method will be those necessary for well-established technical reasons.

# Microbiology of food and animal feeding stuffs — Horizontal method for the enumeration of sulfite-reducing bacteria growing under anaerobic conditions

## 1 Scope

This International Standard specifies a horizontal method for the enumeration of sulfite-reducing bacteria growing under anaerobic conditions. It is applicable to

- products intended for human consumption and the feeding of animals, and
- environmental samples in the area of food production and food handling.

## 2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 6887-1:1999, *Microbiology of food and animal feeding stuffs — Preparation of test samples, initial suspension and decimal dilutions for microbiological examination — Part 1: General rules for the preparation of the initial suspension and decimal dilutions*

ISO 8261, *Milk and milk products — General guidance for the preparation of test samples, initial suspensions and decimal dilutions for microbiological examination*

ISO 7218:1996, *Microbiology of food and animal feeding stuffs — General rules for microbiological examinations*

ISO 7218:1996/Amd.1:2001, *Microbiology of food and animal feeding stuffs — General rules for microbiological examinations — Amendment 1*

ISO/TS 11133-1, *Microbiology of food and animal feeding stuffs — Guidelines on preparation and production of culture media — Part 1: General guidelines on quality assurance for the preparation of culture media in the laboratory*

## 3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

### 3.1

**sulfite-reducing bacteria growing under anaerobic conditions**

bacteria forming countable typical colonies under the conditions specified in this International Standard

## 4 Principle

**4.1** Two agar plates (or tubes) are prepared, using iron sulfite medium, and using a specified quantity of the test sample if the initial product is liquid, or using a specified quantity of an initial suspension in the case of other products.

Two other agar plates (or tubes) are prepared, under the same conditions, using decimal dilutions of the test sample or the initial suspension.

**4.2** The plates (or tubes) are incubated under anaerobic conditions at 48 h 50 °C for to (final reading after 37 °C ± 1 °C 24 h 48 h), or possibly at if thermophilic bacteria are suspected. Typical black-coloured colonies are counted. The black colour of the colonies and the surrounding zone is due to the formation of iron(II) sulfide as a result of the reaction between sulfide ions and trivalent iron [Fe(III)] present in the medium.

**4.3** The number of sulfite-reducing bacteria per millilitre or per gram of sample is calculated from the number of colonies obtained on the plates (or tubes).

**5 Culture medium and diluent**

For current practices, see ISO 7218.

**5.1 Plate count medium: Iron sulfite agar**

<b>5.1.1 Composition</b>		
Enzymatic digest of casein		15 g
Pancreatic digest of soya		5 g
Yeast extract		5 g
Disodium disulfite (Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub> )		1 g
Iron(III) ammonium citrate		1 g
Agar	to a	
	9 g	18 g
Water		1 000 ml

<sup>a</sup> Depending on the gel strength of the agar.

**5.1.2 Preparation**

Dissolve the ingredients in the water by heating.

If necessary, adjust the pH so that after sterilization it is at 7,6 ± 0,2 25 °C.

Pour 250 ml portions of the medium into 500 ml flasks.

If the enumeration is performed by use of tubes (6.5), pour 20 ml or 25 ml of medium into the tubes. Sterilize for 15 min in an autoclave set at 121 °C.

Just before use, de-aerate the medium.

**5.2 Saline peptone diluent**

See ISO 6887-1:1999, 5.2.1.

**6 Apparatus and glassware**

Usual microbiological equipment (see ISO 7218) and, in particular, the following.

**6.1 Homogenization equipment**, for samples of solid food (see ISO 7218).

Two other agar plates (or tubes) are prepared, under the same conditions, using decimal dilutions of the test sample or the initial suspension.

**4.2** The plates (or tubes) are incubated under anaerobic conditions at 48 h 50 °C for to (final reading after 37 °C ± 1 °C 24 h 48 h), or possibly at if thermophilic bacteria are suspected. Typical black-coloured colonies are counted. The black colour of the colonies and the surrounding zone is due to the formation of iron(II) sulfide as a result of the reaction between sulfide ions and trivalent iron [Fe(III)] present in the medium.

**4.3** The number of sulfite-reducing bacteria per millilitre or per gram of sample is calculated from the number of colonies obtained on the plates (or tubes).

**5 Culture medium and diluent**

For current practices, see ISO 7218.

**5.1 Plate count medium: Iron sulfite agar**

<b>5.1.1 Composition</b>		
Enzymatic digest of casein		15 g
Pancreatic digest of soya		5 g
Yeast extract		5 g
Disodium disulfite (Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub> )		1 g
Iron(III) ammonium citrate		1 g
Agar	to a	
	9 g	18 g
Water		1 000 ml

<sup>a</sup> Depending on the gel strength of the agar.

**5.1.2 Preparation**

Dissolve the ingredients in the water by heating.

If necessary, adjust the pH so that after sterilization it is at 7,6 ± 0,2 25 °C.

Pour 250 ml portions of the medium into 500 ml flasks.

If the enumeration is performed by use of tubes (6.5), pour 20 ml or 25 ml of medium into the tubes. Sterilize for 15 min in an autoclave set at 121 °C.

Just before use, de-aerate the medium.

**5.2 Saline peptone diluent**

See ISO 6887-1:1999, 5.2.1.

**6 Apparatus and glassware**

Usual microbiological equipment (see ISO 7218) and, in particular, the following.

**6.1 Homogenization equipment**, for samples of solid food (see ISO 7218).

**6.3 Anaerobic jars**, with equipment for generating an anaerobic atmosphere, and including a system to check the anaerobic conditions.

**6.4 Incubator**, capable of being maintained at  $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$  and, if necessary, at  $50^{\circ}\text{C} \pm 1^{\circ}\text{C}$ .

**6.5 Test tubes**, of dimensions 16 mm × 160 mm, and **flasks** or **bottles** of capacity 500 ml.

## 7 Sampling

It is important that the laboratory receive a sample which is truly representative and has not been damaged or changed during transport or storage.

Sampling is not part of the method specified in this International Standard. If there is no specific International Standard dealing with sampling of the product concerned, it is recommended that the parties concerned come to an agreement on this subject.

## 8 Preparation of test sample

Prepare the test samples in accordance with ISO 6887-1, ISO 8261 or the specific International Standard appropriate to the product concerned. If there is no specific International Standard available, it is recommended that the parties concerned come to an agreement on this subject.

## 9 Procedure

### 9.1 General

A diagram of procedure is given in Annex A.

### 9.2 Test portion, initial suspension and dilutions

See ISO 6887-1, ISO 8261 or any specific International Standard appropriate to the concerned product.

Heat treatment of the initial suspension may be necessary to eliminate vegetative forms of spore-forming bacteria and/or non-spore-forming bacteria. Temperatures and heating times vary according to the actual need, from combinations producing a definite pasteurization effect at a moderate heat activation effect (e.g.  $75^{\circ}\text{C}$  for 20 min), to boiling for several minutes. In this case, results could be given as number of spores of sulfite-reducing bacteria growing under anaerobic conditions.

### 9.3 Inoculation

Take two sterile Petri dishes. Using a sterile pipette, transfer to each dish 1 ml of test sample if the product is liquid, or 1 ml of the initial suspension in the case of other products.

Take two other sterile Petri dishes. Using a fresh sterile pipette, transfer to each dish 1 ml of the first decimal dilution ( $10^{-1}$ ) of the test sample if the product is liquid, or 1 ml of the first decimal dilution of the initial suspension ( $10^{-2}$ ) in the case of other products.

Repeat the described procedure with the further dilutions, using a fresh sterile pipette for each dilution.

Pour into each Petri dish approximately 15 ml of iron sulfite agar (5.1) which has been cooled to  $44^{\circ}\text{C}$  to  $47^{\circ}\text{C}$  in the water bath (6.2). The time elapsing between inoculation of the Petri dishes and addition of the agar should not exceed 15 min. Carefully mix the inoculum with the medium by horizontal movements and allow the medium to solidify.

After the medium has solidified, pour 5 ml to 10 ml of the same medium into the dish as an overlay.



. Mix gently without forming bubbles, and leave the medium to solidify in a cold water bath (6.2).

## 9.4 Incubation


In the case of tubes, incubation in anaerobic jars is not necessary.

The test report shall also state if further tests are to be carried out by a reference laboratory, or, if already carried out, what the results were.




## Annex A (normative)

### Diagram of procedure



## Annex A (normative)

### Diagram of procedure





## BSI — British Standards Institution

BSI is the independent national body responsible for preparing British Standards. It presents the UK view on standards in Europe and at the international level. It is incorporated by Royal Charter.

### Revisions

British Standards are updated by amendment or revision. Users of British Standards should make sure that they possess the latest amendments or editions.

It is the constant aim of BSI to improve the quality of our products and services. We would be grateful if anyone finding an inaccuracy or ambiguity while using this British Standard would inform the Secretary of the technical committee responsible, the identity of which can be found on the inside front cover. Tel: +44 (0)20 8996 9000. Fax: +44 (0)20 8996 7400.

BSI offers members an individual updating service called PLUS which ensures that subscribers automatically receive the latest editions of standards.

### Buying standards

Orders for all BSI, international and foreign standards publications should be addressed to Customer Services. Tel: +44 (0)20 8996 9001.

Fax: +44 (0)20 8996 7001. Email: [orders@bsi-global.com](mailto:orders@bsi-global.com). Standards are also available from the BSI website at <http://www.bsi-global.com>.

In response to orders for international standards, it is BSI policy to supply the BSI implementation of those that have been published as British Standards, unless otherwise requested.

### Information on standards

BSI provides a wide range of information on national, European and international standards through its Library and its Technical Help to Exporters Service. Various BSI electronic information services are also available which give details on all its products and services. Contact the Information Centre. Tel: +44 (0)20 8996 7111. Fax: +44 (0)20 8996 7048. Email: [info@bsi-global.com](mailto:info@bsi-global.com).

Subscribing members of BSI are kept up to date with standards developments and receive substantial discounts on the purchase price of standards. For details of these and other benefits contact Membership Administration. Tel: +44 (0)20 8996 7002. Fax: +44 (0)20 8996 7001. Email: [membership@bsi-global.com](mailto:membership@bsi-global.com).

Information regarding online access to British Standards via British Standards Online can be found at <http://www.bsi-global.com/bsonline>.

Further information about BSI is available on the BSI website at <http://www.bsi-global.com>.

### Copyright

Copyright subsists in all BSI publications. BSI also holds the copyright, in the UK, of the publications of the international standardization bodies. Except as permitted under the Copyright, Designs and Patents Act 1988 no extract may be reproduced, stored in a retrieval system or transmitted in any form or by any means – electronic, photocopying, recording or otherwise – without prior written permission from BSI.

This does not preclude the free use, in the course of implementing the standard, of necessary details such as symbols, and size, type or grade designations. If these details are to be used for any other purpose than implementation then the prior written permission of BSI must be obtained.

Details and advice can be obtained from the Copyright & Licensing Manager. Tel: +44 (0)20 8996 7070. Fax: +44 (0)20 8996 7553. Email: [copyright@bsi-global.com](mailto:copyright@bsi-global.com).